

Human SOAT2/ACACT2/ACAT2 (Trx & His tag) recombinant protein



Catalog Number: 519591

General Information

Gene Name Synonym

Acetyl-CoA acetyltransferase cytosolic, Acetyl-CoA transferase-like protein, Cytosolic acetoacetyl-CoA thiolase, ACTL, ACAT2

Protein Construction

A DNA sequence encoding the Human SOAT2/ACACT2/ACAT2 Met1-Glu397 is expressed with a Trx, 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB□ 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human SOAT2/ACACT2/ACAT2 (Trx & His tag)
recombinant protein