

# GNAQ antibody



Catalog Number: 110991

## Product name

GNAQ antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

GNAQ Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, Jurkat cells, mouse brain tissue, mouse lung tissue, mouse skin tissue. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of GNAQ recombinant protein (Accession Number: NM\_002072). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110991(GNAQ antibody) at dilution of 1:1000