

Catalog Number: 105991

Product name

Anti-NFASC / Neurofascin antibody

Immunogen

[Human NFASC / Neurofascin \(His Tag\) recombinant protein](#)

Specificity

Human NFASC / Neurofascin

Antibody description

Mouse monoclonal to NFASC / Neurofascin

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human NFASC / Neurofascin (rh NFASC / Neurofascin; O94856-12; Met1-Gln939) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of

water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



NFASC / Neurofascin Antibody (PE), Mouse MAb, Flow Cytometry

Flow cytometric analysis of Human NFASC expression on SH-SY5Y cells. Cells were stained with PE-conjugated anti-Human NFASC. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.