

# TAF9B antibody

Catalog Number: 115991

## Product name

TAF9B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TAF9B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive IHC detected in human breast cancer tissue.

## Preparation

This antibody was obtained by immunization of TAF9B recombinant protein (Accession Number: NM\_015975). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IP


## Dilutions

Recommended Dilution:


IP: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115991(TAF9B antibody) at dilution of 1:100



(under 10x lens)

IP Result of anti-TAF9B (IP:Catalog No:115991, 4ug; Detection:Catalog No:115991 1:300) with HeLa cells lysate 4000ug.