

Anti-Annexin V Rabbit antibody



Catalog Number: 166991

Product name

Anti-Annexin V Rabbit antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Annexin V Rabbit polyclonal antibody

Preparation

Antigen: Purified recombinant human Annexin V protein fragments expressed in E.coli.

Formulation

Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/2000

Validations



Western blot detection of Annexin V in 3T3, HeLa, C6, COS7 and MCF7 cell lysates using Annexin V rabbit pAb (1:2000 diluted). Predicted band size: 30 kDa. Observed band size: 30 kDa. Exposure time: 20s.