

Human PLA2G16/HRASLS3/HREV107 (His tag) recombinant protein



Catalog Number: 518991

General Information

Gene Name Synonym

Group XVI Phospholipase A1/A2, Adipose-Specific Phospholipase A2, AdPLA, H-Rev 107 Protein Homolog, HRAS-Like Suppressor 1, HRAS-Like Suppressor 3, HRSL3, HREV107-1, HREV107-3, Renal Carcinoma Antigen NY-REN-65, PLA2G16, HRASLS3, HREV107

Protein Construction

A DNA sequence encoding the Human PLA2G16/HRASLS3/HREV107 Asp12-Asp132 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PLA2G16/HRASLS3/HREV107 (His tag)
recombinant protein