

# DLG3 antibody



Catalog Number: 109991

## Product name

DLG3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DLG3 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 100kd

## Preparation

This antibody was obtained by immunization of DLG3 recombinant protein (Accession Number: NM\_021120). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

IP: 1:1000-1:10000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109991(DLG3 antibody) at dilution of 1:1000

IP Result of anti-DLG3 (IP:Catalog No:109991, 3ug; Detection:Catalog No:109991 1:3000) with mouse brain tissue lysate 4000ug.