

SLC25A17 Polyclonal Antibody



Catalog Number: 162291

Product name

SLC25A17 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SLC25A17

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SLC25A17 (NP_006349.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations

Immunofluorescence - SLC25A17 Polyclonal Antibody

Immunofluorescence analysis of 293T cells using SLC25A17 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.