

SNX3 Polyclonal Antibody



Catalog Number: 164291

Product name

SNX3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SNX3

Preparation

Antigen: A synthetic Peptide of human SNX3

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

Validations



Western blot - SNX3 Polyclonal Antibody

Western blot analysis of extracts of THP-1 cells, using SNX3 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.