

PSAP antibody

Catalog Number: 114291

Product name

PSAP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PSAP Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of PSAP recombinant protein (Accession Number: NM_002778). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114291(PSAP antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:114291(PSAP antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:114291(PSAP antibody) at dilution of 1:100 (under 40x lens)