

# Anti-RPL10L antibody



Catalog Number: 176291

## Product name

Anti-RPL10L antibody

## Specificity

Human, Monkey

## Antibody description

Rabbit polyclonal antibody to RPL10L

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rpl10l. aa range:151-200

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations



Fig1::; Western Blot analysis of various cells using Ribosomal Protein L10L Polyclonal Antibody diluted at 1: 1000