

Anti-ATF3 Rabbit antibody



Catalog Number: 167691

Product name

Anti-ATF3 Rabbit antibody

Specificity

Human

Antibody description

ATF3 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human ATF3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of ATF3 in A549, MCF-7, U251, HL-60 cell lysates using ATF3 Rabbit mAb(1:1000 diluted). Predicted band size: 21kDa. Observed band size: 21kDa.