

# Anti-ATF3 Rabbit antibody



Catalog Number: 167691

## Product name

Anti-ATF3 Rabbit antibody

## Specificity

Human

## Antibody description

ATF3 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human ATF3

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of ATF3 in A549, MCF-7, U251, HL-60 cell lysates using ATF3 Rabbit mAb(1:1000 diluted). Predicted band size: 21kDa. Observed band size: 21kDa.