

Anti-Neurokinin 1 Receptor Rabbit antibody



Catalog Number: 168691

Product name

Anti-Neurokinin 1 Receptor Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Neurokinin 1 Receptor Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Neurokinin 1 Receptor

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/2000-1/10000

IP: 1/20

Validations



Western blot detection of Neurokinin 1 Receptor in CHO-K1,C6,Jurkat cell lysates using Neurokinin 1 Receptor Rabbit mAb(1:500 diluted).Predicted band size:46kDa.Observed band size:46kDa.