

MAGED2 antibody



Catalog Number: 112391

Product name

MAGED2 antibody

IP: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MAGED2 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HepG2 cells, MCF7 cells. Positive IP detected in HepG2 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of MAGED2 recombinant protein (Accession Number: NM_201222). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:112391(MAGED2 antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells, using MAGED2 antibody Catalog No:112391 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

+	-	-	Normal Rabbit IgG
-	+	-	MAGED2 antibody

IP Result of anti-MAGED2 (IP:Catalog No:112391,

MAGED2 antibody



Catalog Number: 112391

4ug; Detection:Catalog No:112391 1:1000) with HepG2 cells lysate 3200ug.