

IL20RB antibody



Catalog Number: 107391

Product name

IL20RB antibody

Specificity

Human; other species not tested.

Antibody description

IL20RB Mouse Monoclonal antibody. Positive WB detected in Hela cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of IL20RB recombinant protein (Accession Number: BC033292). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107391(IL20RB antibody) at dilution of 1:1000