

Anti-SWAP70 (2E9) Mouse antibody



Catalog Number: 168391

Product name

WB, IP

Anti-SWAP70 (2E9) Mouse antibody

Dilutions

Specificity

WB: 1/1000

Human, Mouse, Rat, Monkey

Antibody description

SWAP70 (2E9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human SWAP70 protein fragments expressed in E.coli.

Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG2b

Applications

Validations



Immunoprecipitation analysis of HeLa cell lysate using SWAP70 mouse mAb.



Western blot detection of SWAP70 in COS7,293T,C6,3T3 and HeLa cell lysates and using SWAP70 mouse mAb (1:1000 diluted). Predicted band size: 70kDa. Observed band size: 70kDa.