

Anti-AR antibody



Catalog Number: 141791

Product name

Anti-AR antibody

Specificity

Human Androgen receptor/NR3C4

Antibody description

Mouse Monoclonal to Human AR

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with a synthetic peptide corresponding to the center region of the Human Androgen receptor/NR3C4. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ICC/IF, IF

Dilutions

ICC/IF: 10-50 µg/mL

Validations



Human Androgen receptor/NR3C4
Immunofluorescence(IF) 15375