

# Syntaxin 11 antibody



Catalog Number: 115791

## Product name

Syntaxin 11 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Syntaxin 11 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in human placenta tissue, DU 145 cells, human skeletal muscle tissue, Jurkat cells, L02 cells, mouse lung tissue. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of Syntaxin 11 recombinant protein (Accession Number: XM\_047419437). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:115791(STX11 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:115791(STX11 antibody) at dilution of 1:100 (under 10x lens)