

ADHFE1 antibody



Catalog Number: 107791

Product name

ADHFE1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ADHFE1 Rabbit Polyclonal antibody. Positive WB detected in rat kidney tissue, mouse heart tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 50 kDa, 45 kDa

Preparation

This antibody was obtained by immunization of ADHFE1 recombinant protein (Accession Number: BC064634). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



rat kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107791(ADHFE1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver using Catalog No:107791(ADHFE1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:107791(ADHFE1 antibody) at dilution of 1:50 (under 40x lens)