

C7orf50 antibody



Catalog Number: 108791

Product name

C7orf50 antibody

Specificity

human; other species not tested.

Antibody description

C7orf50 Rabbit Polyclonal antibody. Positive WB detected in L02 cells. Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of C7orf50 recombinant protein (Accession Number: NM_001134396). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:108791(C7orf50 antibody) at dilution of 1:500