

KIF27 antibody



Catalog Number: 112009

Product name

KIF27 antibody

Specificity

Human, Dog; other species not tested.

Antibody description

KIF27 Rabbit Polyclonal antibody. Positive IF detected in MDCK cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM_017014904). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations



Immunofluorescent analysis of MDCK cells using Catalog No:112009(KIF27 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)