

# Anti-SPAK Rabbit antibody



Catalog Number: 168009

## Product name

Anti-SPAK Rabbit antibody

## Specificity

Human

## Antibody description

SPAK Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human SPAK

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC


## Dilutions

WB: 1/1000

IHC: 1/20

FC: 1/20

## Validations



Western blot detection of SPAK in HeLa,A549,HL-60,U2OS cell lysates using SPAK Rabbit mAb(1:1000 diluted).Predicted band size:60kDa.Observed band size:65kDa.