

AMPK beta 2 antibody



Catalog Number: 108009

Product name

AMPK beta 2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

AMPK beta 2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, fetal human brain tissue, HeLa cells, human heart tissue, mouse testis tissue. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of AMPK beta 2 recombinant protein (Accession Number: NM_005399). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108009(PRKAB2 antibody) at dilution of 1:400