

Human SENP7/KIAA1707/SSP2/SUSP2 (His tag) recombinant protein



Catalog Number: 518009

General Information

Gene Name Synonym

Sentrin-Specific Protease 7, SUMO-1-Specific Protease 2, Sentrin/SUMO-Specific Protease SENP7, SENP7, KIAA1707, SSP2, SUSP2

Protein Construction

A DNA sequence encoding the Human SENP7/KIAA1707/SSP2/SUSP2 Met695-Ala864 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM HEPES, 5% Glycerol, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human SENP7/KIAA1707/SSP2/SUSP2 (His tag)
recombinant protein