

Anti-AP-2-alpha Rabbit antibody



Catalog Number: 169009

Product name

Anti-AP-2-alpha Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

AP-2-alpha Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human AP2 alpha

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000-1/5000

ICC/IF: 1/50-1/200

FC: 1/20

IP: 1/20-1/50

Validations



Western blot detection of AP-2-alpha in C6 cell lysates using AP-2-alpha Rabbit mAb(1:1000 diluted).Predicted band size:48kDa.Observed band size:48kDa.