

CD226 antibody

Catalog Number: 109009

Product name

CD226 antibody

Specificity

Human; other species not tested.

Antibody description

CD226 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of CD226 recombinant protein (Accession Number: NM_001303618). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109009(CD226 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109009(CD226 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).