

ERCC3 antibody

Catalog Number: 110409

Product name

ERCC3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ERCC3 Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Positive WB detected in COLO 320 cells, A431 cells, HEK-293 cells, HeLa cells, HT-1080 cells, mouse colon tissue. Observed molecular weight by Western-blot: 89 kDa

Preparation

This antibody was obtained by immunization of ERCC3 recombinant protein (Accession Number: NM_000122). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110409(ERCC3 antibody) at dilution of 1:200



IP Result of anti-ERCC3 (IP:Catalog No:110409, 4ug; Detection:Catalog No:110409 1:500) with A431 cells lysate 1200ug.