

# Nucleoside phosphorylase antibody



Catalog Number: 113409

---

## Product name

Polyclonal

Nucleoside phosphorylase antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

Nucleoside phosphorylase Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, K-562 cells, Raji cells, U-937 cells. Observed molecular weight by Western-blot: 32kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of Nucleoside phosphorylase recombinant protein (Accession Number: NM\_000270). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113409(NP antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality