

RAB11FIP3 antibody



Catalog Number: 114409

Product name

RAB11FIP3 antibody

Immunogen

[Human RAB11FIP3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB11FIP3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 82 kDa

Preparation

This antibody was obtained by immunization of RAB11FIP3 recombinant protein (Accession Number: NM_014700). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114409(RAB11FIP3 Antibody) at dilution of 1:600