

Anti-EPG5 antibody



Catalog Number: 175809

Product name

Anti-EPG5 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to EPG5

Preparation

This antigen of this antibody was recombinant protein of epg5

Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:1000-1:2000

IHC-P: 1:200-1:500

Validations




Fig1:: Western blot analysis of Rat Liver Tissue using EPG5 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at




Fig2:: Immunohistochemical analysis of paraffin-embedded Mouse Liver Tissue using EPG5 Polyclonal Antibody.