

STXBP1 antibody

Catalog Number: 115809

Product name

STXBP1 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

STXBP1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in A375 cells, A549 cells, HeLa cells, human brain tissue, Y79 cells. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 65-68 kDa

Preparation

This antibody was obtained by immunization of STXBP1 recombinant protein (Accession Number: NM_001032221). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:115809(STXBP1 antibody) at dilution of 1:1000



IP Result of anti-STXBP1 (IP:Catalog No:115809, 4ug; Detection:Catalog No:115809 1:1000) with HeLa cells lysate 2400ug.




Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog

STXBP1 antibody

Catalog Number: 115809

No:115809(STXBP1 Antibody) at dilution of 1:50
(under 10x lens)



Immunohistochemistry of paraffin-embedded
mouse brain tissue slide using Catalog
No:115809(STXBP1 Antibody) at dilution of 1:50
(under 40x lens)