

# Anti-Tissue Factor Rabbit antibody



Catalog Number: 167809

## Product name

Anti-Tissue Factor Rabbit antibody

## Specificity

Human

## Antibody description

Tissue Factor Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Tissue Factor

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/1000

IHC: 1/200

## Validations



Western blot detection of Tissue Factor in A549,HL-60,U2OS cell lysates using Tissue Factor Rabbit mAb(1:1000 diluted).Predicted band size:33kDa.Observed band size:45kDa.