

# Cadherin-20 antibody



Catalog Number: 108809

## Product name

Cadherin-20 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Cadherin-20 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human placenta tissue, mouse brain tissue. Observed molecular weight by Western-blot: 120-130kd

## Preparation

This antibody was obtained by immunization of Cadherin-20 recombinant protein (Accession Number: NM\_031891). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108809(CDH20 antibody) at dilution of 1:1000