

General Information

Protein Construction

A DNA sequence encoding the mouse SLAMF8 (NP_083360.2) extracellular domain (Met 1-Asp 231) was fused with the a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to bind biotinylated mouse SLAMF8 in a functional ELISA.

Purity

> 80 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Val 21

Molecular Mass

The secreted recombinant mouse SLAMF8

consists of 222 amino acids and has a predicted molecular mass of 24.8 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rmSLAMF8 is approximately 33-38 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse SLAMF8 / BLAME Protein (His Tag) SDS-PAGE