

NELF-A antibody

Catalog Number: 113109

Product name

NELF-A antibody

IP: 1:200-1:2000

IF: 1:10-1:100


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

NELF-A Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in HeLa cells, Jurkat cells, mouse brain tissue, mouse liver tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 60 kDa




Preparation

This antibody was obtained by immunization of NELF-A recombinant protein (Accession Number: NM_005663). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113109(NELF-A Antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

IP Result of anti-NELF-A (IP:Catalog No:113109, 4ug; Detection:Catalog No:113109 1:500) with mouse brain tissue lysate 3440ug.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IF, WB, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunofluorescent analysis of HeLa cells using Catalog No:113109(NELF-A Antibody) at dilution

NELF-A antibody



Catalog Number: 113109

of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)