

UPF1 Polyclonal Antibody



Catalog Number: 162509

Product name

UPF1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to UPF1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 270-370 of human UPF1 (NP_002902.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - UPF1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using UPF1 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - UPF1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using UPF1 antibody . Blue: DAPI for nuclear staining.