

# MBP-Tag antibody



Catalog Number: 112509

## Product name

MBP-Tag antibody

## Specificity

Tag

## Antibody description

MBP-Tag Rabbit Polyclonal antibody. Positive WB detected in lane1: MBP-NOL6 64kd; lane2 MBP 42kd, Recombinant protein.

## Preparation

This antibody was obtained by immunization of MBP-Tag recombinant protein. Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



lane1: MBP-NOL6 64kd; lane2 MBP 42kd were subjected to SDS PAGE followed by western blot with Catalog No:117331(MBP-Tag Antibody) at dilution of 1:3000