

# ATXN3L antibody



Catalog Number: 108509

## Product name

ATXN3L antibody

## Immunogen

[Human ATXN3L Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ATXN3L Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 41 kDa

## Preparation

This antibody was obtained by immunization of ATXN3L recombinant protein (Accession Number: NM\_001135995). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117070(ATXN3L Antibody) at dilution of 1:600