

Human CTSV/CATL2/CTSL2/CTSU/UNQ268 (His tag) recombinant protein



Catalog Number: 519509

General Information

Gene Name Synonym

Cathepsin L2, Cathepsin U, Cathepsin V, CTSL2, CATL2, CTSU, CTSV

Protein Construction

A DNA sequence encoding the Human CTSV/CATL2/CTSL2/CTSU/UNQ268 Val18-Val334 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM MES, 150mM NaCl, pH 5.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human CTSV/CATL2/CTSL2/CTSU/UNQ268 (His tag) recombinant protein