

# KANK2 antibody



Catalog Number: 111909

## Product name

KANK2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KANK2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse lung tissue, mouse heart tissue, mouse kidney tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 95kd

## Preparation

This antibody was obtained by immunization of KANK2 recombinant protein (Accession Number: NM\_001379563). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111909(KANK2 antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using KANK2 antibody Catalog No:111909 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).