

# MVB12A Polyclonal Antibody



Catalog Number: 162909

## Product name

MVB12A Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to MVB12A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human MVB12A (NP\_612410.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - MVB12A Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using MVB12A antibody at 1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.