

# RPS4Y1 antibody



Catalog Number: 114909

## Product name

RPS4Y1 antibody

## Immunogen

[Human RPS4Y1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RPS4Y1 Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells, A2780 cells, HeLa cells, mouse ovary tissue, mouse testis tissue, SKOV-3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 29-32kd

## Preparation

This antibody was obtained by immunization of RPS4Y1 recombinant protein (Accession Number: NM\_001008). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:114909(RPS4Y1 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using RPS4Y1 antibody Catalog No:114909 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).