

Anti-MERTK Rabbit antibody



Catalog Number: 167909

Product name

Anti-MERTK Rabbit antibody

Specificity

Human

Antibody description

MERTK Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MERTK

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/50-1/500

IP: 1/20

Validations



Western blot detection of MERTK in MCF-7,293T cell lysates using MERTK Rabbit mAb(1:1000 diluted).Predicted band size:110kDa.Observed band size:210kDa.