

AFM antibody



Catalog Number: 107909

Product name

AFM antibody

Specificity

Human; other species not tested.

Antibody description

AFM Rabbit Polyclonal antibody. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 85 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001133). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:107909(AFMM antibody) at dilution of 1:1000