

# Human UBE2D4/UBCH5D (GST tag) recombinant protein



Catalog Number: 518909

## General Information

### Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 D4, HBUCE1, Ubiquitin Carrier Protein D4, Ubiquitin-Protein Ligase D4, UBE2D4, UBCH5D

### Protein Construction

A DNA sequence encoding the Human UBE2D4/UBCH5D Met1-Met147 is expressed with a GST tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human UBE2D4/UBCH5D (GST tag) recombinant protein