

# Human MCL1/BCL2L3 (His tag) recombinant protein



Catalog Number: 519909

## General Information

### Gene Name Synonym

Induced myeloid leukemia cell differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl2-L-3; Bcl-2-related protein EAT/mcl1; mcl1/EAT; MCL1; BCL2L3

### Protein Construction

A DNA sequence encoding the Human MCL1/BCL2L3 Met1-Gly327 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 200mM NaCl, 2mM DTT, 20% glycerol, pH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human MCL1/BCL2L3 (His tag) recombinant protein