

Anti-OR5A2 antibody

Catalog Number: 176209

Product name

Anti-OR5A2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR5A2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or5a2. aa range:267-316

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of 3T3 cells using Olfactory receptor 5A2 Polyclonal Antibody