

# Anti-SMARCC1/BAF155 (2A4) Mouse antibody



Catalog Number: 168209

## Product name

Monoclonal

Anti-SMARCC1/BAF155 (2A4) Mouse antibody

## Ig Type

## Specificity

Mouse IgG1

Human

## Applications

## Antibody description

WB

SMARCC1/BAF155 (2A4) Mouse monoclonal antibody

## Dilutions

WB: 1/1000


## Preparation

Antigen: Purified recombinant human SMARCC1 protein fragments expressed in E.coli.

## Validations

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of SMARCC1 in K562, Jurkat and HeLa cell lysates using SMARCC1 mouse mAb (1:1000 diluted). Predicted band size: 155kDa. Observed band size: 155kDa.

## Clonality