

Human LOXL2 (Avi & His tag) recombinant protein



Catalog Number: 518209

General Information

Gene Name Synonym

Lysyl oxidase homolog 2; Lysyl oxidase-like protein 2; Lysyl oxidase-related protein 2; Lysyl oxidase-related protein WS9-14; LOXL2

Protein Construction

A DNA sequence encoding the Human LOXL2 Gln26-Gln774 is expressed fused with a Avi-His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 50mM Sodium borate, 10mM CaCl₂, 1.2M Urea, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human LOXL2 (Avi & His tag) recombinant protein